

ClearKlens Cleansinald

General purposes surface detergent and disinfectant

Description

ClearKlens Cleansinald is a concentrated detergent disinfectant for the pharmaceutical and cosmetic industries.

ClearKlens Cleansinald is a surface active product based on alkylamines and quaternary ammonium.

Key properties

- **ClearKlens Cleansinald** is used for cleaning and disinfection of walls, floors, work surfaces and process equipment in pharmaceuticals and cosmetics industries, and for animal research centres.
- The formula of **ClearKlens Cleansinald** is a combination of quaternary ammonium compounds and alkylamines.
- **ClearKlens Cleansinald** has a broad spectrum of activity.
- **ClearKlens Cleansinald** can be applied as a spray, wiping or wet sweeping without rinse.
- Small items or dismantled equipment may be immersed in a **ClearKlens Cleansinald** solution.

Benefits

- No aldehyde
- Bactericidal, yeastcidal and virucidal activity
- Efficient on fat soils and proteins
- No rinse required
- No effect on materials normally encountered

Use instructions

ClearKlens Cleansinald is usually used at 0.5% v/v aqueous solution (20mL for 4L of cold or lukewarm water).

ClearKlens Cleansinald is used without rinse.

VH9



Pharma™

ClearKlens Cleansinald

General purposes surface detergent and disinfectant

Technical data

Appearance	Clear, slightly yellow liquid
pH (1% 20°C)	11
Density (20°C)	1.05
Chemical Oxygen Demand (COD)	563 gO ₂ /kg
Nitrogen content (N)	14 g/kg
Phosphorus content (P)	None

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Store in original closed containers or (where applicable) in approved bulk tanks, away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is provided in a separate Material Safety Data Sheet.

Product compatibility

Applications solutions, when used as directed, will not affect materials normally encountered in the aforementioned industries.

Test method

Each lot is provided with:

Certificate of Analysis

Environmental information

Registration

Microbiological data

EN 1276 bactericidal: passed at 0.75% dilution against *E. coli*, *P. aeruginosa*, *S. aureus* and *E. hirae* in hard water (300ppm as CaCO₃), high soil (0.3% bovine albumin) and 5minutes contact time at 20°C.

EN 1650 yeasticidal: passed at 0.5% against *C. albicans* in hard water (300ppm as CaCO₃), in high soil conditions (0.3% bovine albumin respectively) with a 15 minute contact time at 20°C.

EN 13697: passes passed at 0.5% dilution in hard water (300ppm as CaCO₃) in low soil conditions (0.03% bovine albumin) with a 5 minute contact time and passes against *C. albicans* in hard water (300ppm as CaCO₃), in low soil conditions (0.03% bovine albumin) with a 15 minute contact time.

ClearKlens Cleansinald possesses virucidal activity in no soil conditions at 20°C against Hepatitis B and HIV after 5 minutes contact time at 0.25%.

ClearKlens Cleansinald possesses virucidal activity in no soil conditions at 20°C against Rotavirus after 15 minutes contact time at 0.5%.

Full microbiological reports available upon request.

Article code	Pack size	Language
7513417	2 x 5L	UK, IE, FR, AT, SE, FI, GR, SP